



Disposal Sets Overview

Sets for the safe disposal of liquid and explosive laboratory waste.

- ✓ Electrically conductive material to avoid explosion hazards.
- ✓ Sets with a practical lid or solid ball valve funnel and grounding cables.
- ✓ Safety canister with mechanical level control.
- ✓ Optionally with collecting tray to protect against leaks.
- ✓ All sets can be used immediately without additional accessories.

All set components can also be reordered individually:



Component	Funnel			Canister			Collection container	
	Type	Model	Feature	Content	Thread	Part number	Set	
Type	Ball valve	Ball valve	Hinged lid	Safety canister	Safety canister	Safety canister	Collecting tray	
Model	ARNOLD	ARNOLD	Marco	-	-	-	-	
Feature	Lance (rem.)	with lance (220 mm)	without lance	mechanical level control	-	mechanical level control	with base insert	
Content	-	-	-	10 litre	10 litre	20 litre	10 litre	
Thread	S60/61	S60/61	S60/61	S60/61	S60/61	S60/61	-	
Set	Part number	317 621	317 621	318 960	108 042	107 953	108 043	117 985
A	307 320	•			•			
B	307 315			•		•		•
C	307 416			•	•			•
D	307 321		•				•	